

1 mL  
Syringe  
now  
available



# eVol<sup>®</sup> XR - Every One an Expert

- Speed up and simplify laboratory workflow
- Improve accuracy and reproducibility
- Standardize results independent of operator skill



 **SGE Analytical Science**  
[www.sge.com](http://www.sge.com)

AUSTRALIA & PACIFIC REGION

CHINA

EUROPE

JAPAN

USA

## eVol® is the world's first digital analytical syringe

eVol® is the coupling of two precision devices: a digitally controlled electronic drive and an XCHANGE® enabled analytical syringe.

- XCHANGE® syringes are easily and quickly changed allowing them to be dedicated to individual liquids or methods to prevent possible cross-contamination of reagents.
- eVol® is ergonomic, comfortable and easy to use.
- eVol® is easily calibrated and calibration factors saved for each syringe, enabling laboratories to comply with stringent global laboratory standards (e.g. GLP, GMP, FDA).
- eVol® is suitable for use with volatile samples.
- eVol® is programmable to store a laboratory workflow (up to 98 steps).
- eVol® is suitable for direct injection onto a chromatography column with a consistent flow rate.
- eVol®'s stainless steel needle enables direct injection through septa.
- Password protection options enabling standardization of work processes
- Variable speed of operation


## eVol® applications

Applications for eVol® include:

- Preparation of calibration standards
- Preparation and addition of internal standards
- Precise dispensing of aqueous and non-aqueous liquids
- Routine dispensing
- Ergonomic operation with substances in a fume hood
- Sample dilution
- Instrument (GC and LC) injections
- Eliminating the need for Serial dilutions
- Micro titrations
- TLC spotting
- FDA methods requiring a 1 mL syringe

New  
1 mL Syringe  
expands  
applications  
for eVol®

## Volume range of eVol® XCHANGE® syringes



The award winning eVol® improves the pace of laboratory processes while delivering improved accuracy and reproducibility.

The table below includes examples of how eVol® improves standard laboratory processes

| PROCESS                           | Without eVol®   | With eVol®   | eVol® BENEFIT  |
|-----------------------------------|---|--|--|
| Standard Preparation              | Standards prepared in a large volume flask. From this standard aliquots are individually dispensed into autosampler vials.  | Standards are made up directly in the vial, including the make up solvent.   | <ul style="list-style-type: none"> <li>• Less glassware usage</li> <li>• Reduces waste fluid</li> <li>• Significant time saving</li> <li>• Improved accuracy and reproducibility</li> </ul>  |
| Addition of Standards             | Small amounts of standard aspirated and dispensed into all samples before being transferred to an autosampler vial.   | One aspiration and a fast series of repeated accurate dispenses directly into vials.   | <ul style="list-style-type: none"> <li>• Significant time saving</li> <li>• Improved accuracy and reproducibility</li> </ul>   |
| Delivery of derivatization agents | Laboratory staff required to work in a fume hood with potentially hazardous materials, to prepare combinations of derivatization agents in open vials.  | Process completed with eVol® programmed to aspirate an amount of solvent or agent and then dispense aliquots into sealed vials, Single handed operation. | <ul style="list-style-type: none"> <li>• Improved operator safety, lower spill and splash risk</li> <li>• Ergonomic benefits behind fume hood screen</li> <li>• Improved accuracy and reproducibility</li> <li>• Less glassware use</li> </ul> |
| Serial dilutions                  | Transfer of a small amount of solution to another container. Solvent added to achieve the required volume. This is repeated multiple times to obtain the required final accurate concentration. | One aspiration of the solution can be dispensed directly into the solvent to achieve the required accurate concentration.                                | <ul style="list-style-type: none"> <li>• Complete workflow simplification</li> <li>• Significant time savings</li> <li>• Improved accuracy</li> <li>• Less solvent required</li> <li>• Less glassware used</li> </ul>                          |

### eVol® NMR Edition

eVol® NMR Edition features long stainless steel needles enable 'in tube' sample dilution and mixing, and recovery of samples facilitating re-use of NMR tubes. Contact SGE for details.

### eVol® MEPS™

The eVol® custom programming function enables MEPS™ (micro SPE) to be semi-automated and proofed before transition to fully automated platforms. Ask SGE for more information on how eVol® MEPS™ can improve your method development.



## eVol® XR Digitally Controlled Syringe - Starter Kit

| Description  | Part No. |
|--|----------|
| <b>eVol® XR Syringe Starter Kit</b><br><br>Contains: <ul style="list-style-type: none"> <li>eVol® XR Electronic Syringe</li> <li>3 eVol® Syringes – 5 µL, 100 µL and 1 mL</li> <li>Stand</li> <li>Universal Charger</li> <li>Comprehensive Instruction Manual</li> <li>Disc with Manual in Multiple Languages</li> </ul>   | 2910200  |
| <b>eVol® NMR Edition</b><br><br>Contains: <ul style="list-style-type: none"> <li>eVol® Electronic Syringe</li> <li>3 eVol® Syringes – 5 µL syringe is supplied with a 115 mm needle, 50 µL and 500 µL syringes are supplied with both 115 and 180 mm needles</li> <li>Stand.</li> <li>Universal Charger</li> <li>Comprehensive Instruction Manual</li> <li>Disc with Manual in Multiple Languages</li> </ul> | 2910100  |
| eVol® Electronic Syringe (handle only)   | 2910005  |

## eVol® XCHANGE® and eVol® MEPS™ Syringes

| Description                      | Needle Length (mm) | Needle Gauge | Needle OD (mm) | Needle ID (mm) | Needle Tip | Replacement Needle Part No. | Replacement Plunger Part No. | Syringe 3 Pack Part No. | Syringe Part No. |
|----------------------------------|--------------------|--------------|----------------|----------------|------------|-----------------------------|------------------------------|-------------------------|------------------|
| 5 µL                             | 50                 | 25           | 0.5            | 0.12           | Bevel      | 036910                      | 2910380                      | 2910320                 | 2910020          |
| 5 µL (supplied without needle)   | –                  | –            | –              | –              | –          | –                           | 2910380                      | -                       | 2910021          |
| 50 µL                            | 50                 | 25           | 0.5            | 0.2            | Bevel      | 038110 <sup>Y</sup>         | 2910382                      | 2910322                 | 2910022          |
| 50 µL (supplied without needle)  | –                  | –            | –              | –              | –          | –                           | 2910382                      | -                       | 2910023          |
| 50 µL for MEPS™ applications*    | –                  | –            | –              | –              | –          | –                           | 2910382                      | -                       | 2910027          |
| 100 µL                           | 50                 | 25           | 0.5            | 0.2            | Bevel      | 038110                      | 2910383                      | 2910329                 | 2910029          |
| 100 µL for MEPS™ applications*   | –                  | –            | –              | –              | –          | –                           | 2910383                      | –                       | 2910028          |
| 500 µL                           | 50                 | 23           | 0.63           | 0.32           | Bevel      | 039110 <sup>A</sup>         | 2910384                      | 2910324                 | 2910024          |
| 500 µL (supplied without needle) | –                  | –            | –              | –              | –          | –                           | 2910384                      | -                       | 2910025          |
| 500 µL for MEPS™ applications*   | –                  | –            | –              | –              | –          | –                           | 2910384                      | -                       | 2910026          |
| 1 mL <sup>#</sup>                | 50                 | 23           | 0.63           | 0.32           | Bevel      | 039110                      | 2910385                      | 2910335                 | 2910035          |

<sup>#</sup> 1 mL syringe can only be used with eVol® XR, not compatible with eVol® Classic. <sup>Y</sup> All SGE 25 – 500 µL replacement needles can be used with 50 µL and 100 µL eVol® syringes. <sup>A</sup> All SGE 1 – 2.5 mL replacement needles can be used with 500 µL eVol® syringes. \* The eVol® MEPS™ syringes can be used with the range of MEPS™ BINs.

### Accessories for eVol® XR



Part No. 2910010  
eVol® Stand



Part No. 2910012  
eVol® Universal Charger



Part No. 2910040  
eVol® Replacement Battery



Part No. 2910030  
eVol® Charging Stand

#### AUSTRALIA & PACIFIC REGION

SGE Analytical Science Pty Ltd  
Toll Free: 1800 800 167  
Tel: +61 (0) 3 9837 4200  
Fax: +61 (0) 3 9874 5672  
Email: support@sge.com

#### CHINA

SGE Shanghai Representative Office  
Tel: +86 21 6407 9382  
Fax: +86 21 6407 9386  
Email: china@sge.com

#### EUROPE

SGE Europe Ltd  
European Head Office  
Toll Free: 00800 2790 8999  
Toll Free Fax: 00800 2626 2609  
Tel: +44 (0) 1908 568 844  
Fax: +44 (0) 1908 566 790  
Tel France: +33 (0)1 69 29 80 90  
Fax France: +33 (0)1 69 29 09 25  
Tel Germany: +49 (0) 6155 / 60746 0  
Fax Germany: +49 (0) 6155 / 60746 50  
Email: europe@sge.com

#### JAPAN

SGE Japan Inc  
Tel: +81 45 222 2885  
Fax: +81 45 222 2887  
Email: japan@sge.com

#### UNITED STATES OF AMERICA

SGE Incorporated  
Toll Free: (800) 945 6154  
Tel: +1 512 837 7190  
Fax: +1 512 836 9159  
Email: usa@sge.com

